

Searching for Improved Treatment Options for Human Internal Decontamination

A Nanoscale Approach

Michael D. Kaminski,^{a,c} Carol J. Mertz,^c Martha R. Finck,^c Yumei Xie,^d Haitao Chen,^d Vivian Sullivan,^c Sandra Guy,^b Vincent Turitto,^d and Axel J. Rosengart,^{a,b}
^aCo-Principal Investigator, ^bPritzker School of Medicine, The University of Chicago, ^cArgonne National Laboratory, ^dIllinois Institute of Technology

Overall Research Goal
Develop, design, and demonstrate a novel, integrated system based on magnetic nanoparticles for selective, rapid removal of biological, chemical, and radioactive biohazards from humans.

Current Detoxification Methods and Limitations
Hemodialysis and Hemofiltration: Long procedure duration, extracorporeal circulation of large blood volumes, large-bore arterial access, non-selective substance removal, limited to hydrophilic substances with lower molecular weight
Plasmapheresis: Generally restricted to autoimmune diseases
Extracorporeal Immunoabsorption: Specific removal method but less effective, restricted to specific antibody-antigen interactions, requires circulation of large blood volumes
Direct Injection of Chelators and Antibodies ("Immunotherapy"): Incomplete antigen binding, relatively high antibody dosing; bound toxin can still be toxic or dissociate, leading to rebound intoxication


Overview of Proposed System

Injection
Magnetic nanospheres are developed that are biostabilized with poly(ethylene glycol) (PEG) to avoid rapid biodegradation. They are biodegradable and non-toxic and present a variety of specific toxin-binding receptors. These nanospheres are injected directly into the blood stream of exposed humans.



Chelation
Once injected, the nanospheres circulate freely through the blood stream, offering selective capture and sequestration of blood-borne toxins to the surface receptors.


Separation
The magnetic separator device allows blood to pass through a closed-loop system containing tubular micro-channels designed to optimize blood flow and particle removal. Small magnets within the device trap the particle-toxin complexes at the tube walls and the detoxified blood is returned into the body. After the procedure, the toxin-bound particles can be flushed from the removable tubing and used as an analyte in further testing.

Advantages
Therapeutically Superior: Active removal of biotoxins from exposed humans, not merely binding and passive secretion
Biocompatible: Non-toxic, biodegradable nanoparticles not removed will pose no harm to the subject.
Diverse: Already existing and newly designed antitoxins, antibodies, and ligands can be attached and multiple toxins can be removed ('cocktail' approach)
Repeatable: Chronic exposures or exposures with high tissue levels can be retreated
Compact: Total MNP-antitoxin injected expected to be <2 mg/kg body weight; injectants and magnetic separator device engineered as hand-held, single-use, pre-sterilized, self-applicable unit (in-field use, etc.)
Concentrates Toxin: Collected toxin provides a highly concentrated analyte for further identification.





Magnetic Nanoparticles

**Controlling Size**
Oil-in-water emulsion (copolymer-magnetite-PLA-chloroform) has proven to be a poor technique. Membrane emulsification promises an improvement. The surfactant properties of PEG reduce emulsion size from 1-4 μm to 275 nm


**Vascular Survival of PEGylated Nanospheres**
Poloxamers
• PEG-PPG-PEG triblock polymers
• PPG: poly(propylene glycol) (hydrophobic)
• PEG: poly(ethylene glycol) (hydrophilic)


Attaching Receptors to PEG Terminus
The literature describes a coblock polymer synthesis with biotin end groups. These provide a direct method of attaching streptavidin or other receptors to the particle. This polymer is now being synthesized in our laboratory.


Surface Charge
A charged particle surface attracts plasma proteins; hence the surface charge needs to be near neutral (at pH 7). PEG neutralizes the surface and its flexible chains reduce the immune response.



Magnetic Separation

**High magnetization particles**

**Magnetic field gradients**

**In Vitro Separation of HRP**

In vitro tests were conducted in both 0.9% saline and blood using 10 mg functionalized nanospheres per 40 mL fluid. The tests showed rapid kinetics and efficient removal (~70% in 15 min) of surrogate toxin (HRP) using magnetic nanoparticles.

Contact Time	15 min	60 min
Injected HRP (mg/mL)	0.38	0.38
Bound HRP (mg/mL)	0.27	0.27
Free HRP (mg/mL)	0.11	0.10
% HRP Bound	71.9	72.8