

Laser Measurements for High-Pressure Combustion: Challenges and Opportunities

Robert P. Lucht

School of Mechanical Engineering

Purdue University

West Lafayette, Indiana

Presentation at the Workshop on Techniques for High-Pressure Combustion

Argonne National Laboratories

August 30, 2011


Acknowledgments

- For dual-pump CARS, funding from NASA Glenn under Cooperative Agreement Number NNX07AC90A , technical discussions with Drs. Yolanda Hicks, Clarence Chang, and Randy Locke. For fs CARS, funding from DOE Office of Basic Energy Sciences, AFOSR and NSF.
- Graduate students: Mathew Thariyan, Ning Chai, Daniel Richardson, Aizaz Bhuiyan
- Research staff: Scott Meyer, Yu Matsutomi, Sameer Naik
- Colleagues: Profs. Hukam Mongia, Jay Gore (Purdue), Sukesh Roy and Waruna Kulatilaka (Spectral Energies), Jim Gord (AFRL)

Outline of the Presentation


- **Optically Accessible Gas Turbine Combustor Facilities**
- **Dual-Pump CARS Measurements: Challenges and Optical System**
- **Femtosecond CARS for Single-Shot Temperature at 5-10 kHz Data Rates**
- **Laser System for 5-10 kHz PLIF and PIV, 5 kHz OH PLIF Measurements**
- **Conclusions**

Purdue Gas Turbine Combustion Facility (GTCTF)




Purdue Gas Turbine Combustion Facility (GTCF)

High Pressure Lab System	Maximum Flow Capacity	Max Operating Condition
Natural Gas Heated High Pressure Air	9 lbm/sec	700 psi / 540 deg C 1000 deg F
Electric Heated Air or Nitrogen	1 lbm/sec	600 psi / 600 deg C
Nitrogen	2 to 5 lbm/sec	1,500 psi
Liquid Aviation Fuel (Kerosene)	1 lbm/sec/tank (2 tanks)	1,500 psi
Natural Gas	1 lbm/sec	3500 psi




NASA 9-Point LDI Assembly (Top-Hat)

- Nine simplex injectors arranged at throats of nine converging-diverging venturis in a 3 x 3 arrangement.
- Axial swirlers with helical vanes at 60° impart swirl to incoming heated air.
- Only central injector used for testing.





Coherent Anti-Stokes Raman Scattering (CARS)


- Conventional “Single-Pump” CARS
- Noninvasive
- Coherent Laser-Like Signal
- Spatially and Temporally Resolved
- Excellent Gas Temperature Data (especially at higher temperatures)



Dual-Pump CARS of N₂/CO₂




Overall Experimental System




Measurement Challenges in GTCF

- Translation of probe volume inside the flame zone.
- Installation of pin-hole for spatial overlap of CARS beams not possible, must use reference leg alignment.
- Measurement of non-resonant signal in the reference leg for spectral normalization of CARS signal.
- Safety of thin window, CARS beams are focused tightly in the middle of the test section.

- A:** Aperture
BD: Beam Dump
CL: Camera Lens
IF: Interference Filter
L: Lens
1/2: Half-wave-plate
Pol: Polarizer
S: Shutter
TS: Translation Stage
TL: Zoom Lens



Optical System near GTCF





Operating Conditions, Measurement Locations and Sample DP-CARS Spectra

	$\Phi=0.4$	$\Phi=0.59$	$\Phi=0.80$	$\Phi=1.0$
100 psia (7.0 atm.)	■	■	■	■
125 psia (8.5 atm.)	■			
150 psia (10 atm.)	■			


- Burner Inlet
Temperature: 850 °F (725 K)
- Fuel: Jet-A
- Normalized injector pressure drop = 4%

$F = 0.6$
 $P = 100$ psi
 $DP/P = 4\%$


Note: Distance between points along the centerline is 5 mm



Measurement Grid for DP-CARS




Flame Characteristics @ 100 psia



$\phi = 0.4$



$\phi = 0.59$







$\phi = 0.8$



$\phi = 1.0$

Temp PDFs Along Centerline


Combustor Pressure: 104 psia, Equivalence Ratio: 0.4



Temperature and CO₂/N₂ Correlation Plot

Combustor Pressure: 150 psia., Equivalence Ratio: 0.48

Temp. vs CO₂/N₂ Correlation: Z= 30 mm, Y= 0 mm



----- Adiabatic Equilibrium

Accomplishments and Conclusions

- A new OPO/PDA system was used to generate the 560-nm pump beam in the dual-pump CARS system. Considerable care in laser system alignment was required to obtain good beam quality in the combustor test cell.
- The Zaber translation stages performed well, alignment was maintained over the entire spatial region of interest during the test.
- The reference leg was invaluable for alignment and for frequent recording of the nonresonant signal. Alignment was maintained before and after translation of the large 2-inch prisms.


Accomplishments and Conclusions

- Estimated uncertainty in temperature measurements :
 - Accuracy: 1-2%
 - Precision: 2-3%
- Uncertainty in CO₂/N₂ ratio measurements :
 - Very dependent on CO₂ concentration and on the temperature, approximately 10% relative standard deviation in the range of 5% CO₂ concentration around 1500 K.
- Probe volume dimensions:
 - 500-600 μm (FWHM) along the laser propagation direction.
 - 50 μm perpendicular to the laser direction.

Modified Combustor Window Assembly

- Cross section increased from 3"x3" to 4.2"x4.2". The modified CWA is fabricated from Hastelloy-X instead of stainless steel. Brazing has been eliminated. Film cooling air passages are incorporated in the injector assembly rather than in the CWA. Thermal barrier coatings are being applied to the window assembly inner surfaces.
- Upstream spool section has been redesigned to accommodate the larger injectors and to ensure uniform flow into the injector.
- Downstream spool sections redesigned for larger flow cross section.


Test Rig with Modified Combustor Window Assembly (CWA) Installed




Single-Shot CPP fs CARS

Nonresonant spike



+1 ps Probe Delay



20 Single-Shots with Chirped Probe Pulse




Temperature Histograms from Single-Shot fs CARS in Flames



Single-Shot Temperature Measurements at 5 kHz

- New laser system from Coherent delivered to Purdue (AFOSR DURIP Program).
- Laser rep rate of either 5 kHz with 2.6 mJ per pulse or 10 kHz with 1.0 mJ per pulse.
- Pulse durations of either 60 fs or 30 fs.
- Pulse shaper integrated into the system between the oscillator and amplifier.


Single-Shot Temperature Measurements at 5 kHz: Coherent Ultrafast Laser System





Single-Shot Temperature Measurements at 5 kHz

- Demo camera systems loaned by Princeton Instruments and Andor.
- Princeton Instruments system: 512x512 chip, physical mask, spectral acquisition from 66 rows by 512 columns, 5 kHz or 10 kHz spectral acquisition with charge in columns vertically binned.
- Andor system: 1024x1024 chip, optical mask, 50 rows by 1024 columns, 5 kHz spectral acquisition with charge in columns vertically binned, also have 512x512, 10 kHz system.

Single-Shot Temperature Measurements at 5 kHz







Single-Shot Temperature Measurements at 5 kHz







Single-Shot Temperature Measurements at 5 kHz


Hencken Burner Calibration Flames



Single-Shot Temperature Meas at 5 kHz: Nonresonant Background Suppressed



Single-Shot Temperature Meas at 5 kHz: Nonresonant Background Suppressed



Conclusions: Fs CARS

- Temperature determined from single-shot fs CARS N₂ spectra recorded at 1 and 5 kHz from laminar, unsteady, and turbulent flames.
- Signals are very strong and reproducible from shot to shot. Precision ~~comparable to or~~ better than the best single-shot ns CARS systems.
- Theoretical model developed to fit CPP fs CARS spectra for temperature, accuracy of measurements can be improved. Main issue is incomplete characterization of the laser beams.
- CPP fs CARS concentration measurements in progress. Data analysis issues similar as for measurement of temperature.


Future Work: fs CARS

- Temperature measurements in turbulent flames.
- Comparison of measurements with and without background suppression.
- Concentration measurements are still an open issue.
- More accurate characterization of the ultrafast laser pulses.
- Measurements in high pressure flames.
- Electronic resonance enhancement.

High-Repetition-Rate Diagnostic Techniques Based on Diode-Pumped Nd:YAG Lasers

- 10 kHz PIV – dual-head Edgewave laser
- 10 kHz OH PLIF – Credo dye laser pumped by the Edgewave laser


High-Repetition-Rate Laser System



Edgewave Diode-Pumped Solid State Nd:YAG Laser: 5 kHz Rep Rate, Dual-Head, 6 mJ/Pulse at 532 nm, 7 nsec Pulses


Sirah Credo Dye Laser
5 kHz Rep Rate, 500 mJ/
Pulse at 283 nm (2.5 W
average power in UV)

5 kHz OH PLIF of Liquid and Gelled Droplet Combustion




- Variety of gelled and liquid droplets have been investigated – methanol and JP8 fuels, HPC and aerosil gellants
- Drop sizes on the order of a few mm
- Combustion process is much more dynamic for the gelled droplets
- Effects of pressure investigated in optical cell


5 kHz OH PLIF of Liquid and Gelled Droplet Combustion



5 kHz OH PLIF of Liquid and Gelled Droplet Combustion



5 kHz OH PLIF of Liquid and Gelled Droplet Combustion



High-Pressure Diagnostic Techniques

- **Measurements in high-pressure flames complicated because of the experimental apparatus (windows, etc.) and because of increasing collisional rates (increased quenching, absorption line widths). High-rep-rate PLIF and PIV will be extremely valuable for turbulent flames.**
- **CARS has been demonstrated in high-pressure flames systems for temperature and concentration. Fs CARS offer the potential for high-rep-rate single-shot measurements in high-pressure systems.**
- **Dual-resonance techniques like electronic-resonance-enhanced CARS or two-color polarization spectroscopy may help to overcome interferences due to overlapping absorption lines from different species.**